

# Magister C24 with Cellbind



	Column agglutination techniques	Immunofixation with column technology
Easy to interpret	✓	✓ ✓
Reagents free of position	—	✓ ✓
Clear and sharp results	✓	✓ ✓
Minimal risk on false interpretation	—	✓ ✓
Allowing optimal digital processing of data	✓	✓ ✓
Minimal occurrence of smear in column	—	✓ ✓



#### References

1. Wan, H-L et al. Poster presentation ISBT; 2012
2. Lee, C-C et al. Poster presentation ISBT; 2012
3. Donk, E.M.M. van der, et al. Internal and external trial; 2013
4. Waldekker, F. Thesis; 2010
5. Sanquin, data on file
6. Cellbind Screen and Magister in the daily routine of the Clinical Laboratory: Poster ISBT 2011, N. Osmanovic, et al



**Magister C24 & Cellbind,  
the ideal combination**



# MORE FLEXIBILITY

#### Sanquin Reagents BV

Plesmanlaan 125  
NL-1066 CX Amsterdam  
P.O. Box 9190  
NL-1006 AD Amsterdam

[www.sanquin.org/reagents](http://www.sanquin.org/reagents)

January 2018  
K-0001-01

#### For order enquiries

Telephone +31 20 - 512 37 39  
Fax + 31 20 - 512 35 70  
E-mail [reagentsorders@sanquin.nl](mailto:reagentsorders@sanquin.nl)

#### For technical information

Telephone +31 20 - 512 35 99  
E-mail [reagents@sanquin.nl](mailto:reagents@sanquin.nl)

#### Sanquin:

Blood Bank  
Diagnostic Services BV  
Plasma Products BV  
Research  
Reagents BV

For more information, pls contact:

Tel nos: 441-0399 / 441-0400

Email: [marketing@mrlcybertec.com](mailto:marketing@mrlcybertec.com)

# MAGISTER C24

## first compact fully automated blood group analyser

Magister C24 ensures safe, accurate and efficient blood group analysis

### Magister C24 is flexible and fast:

- Fully automated
- Table top, low footprint
- Continuous loading
- Easy to load
- Stat samples
- Connected to LIS
- Modular, hence low maintenance time



**FLEXIBLE  
ECONOMIC  
SAFE**

# CELLBIND

## immunofixation column technology

One type of card for all tests



**FIT  
FOR THE  
FUTURE!**

### Cellbind Screen:

- Limited stock and storage space required
- Next sample on the same card: use all columns

### Cellbind Screen suitable for:


- Screening and antibody identification
- Typing of blood group antigens
- Reverse typing
- Crossmatch
- Direct Antiglobulin Test
- Antibody titration

### Magister C24 and Cellbind



a single card for any protocol > minimal stock required  
card is easy to finetune to your protocol > suits your own designed protocol  
results sharp readable above and below in gel > precise and clear results

### benefits for your practice:


software interprets reactions of the columns > human mistakes are prevented  
next sample may start on the very next available column > no spillage  
protocols are easy to adapt > ideal when diagnostic guidelines are altered



**Optional incubation 15 minutes 37°C**



**Centrifugate for 10 minutes**



**positive      negative**

