

EMPOWERING PHILIPPINE LABORATORIES with MUREX DIAGNOSTICS SOLUTIONS

Hospital. Clinic. Emergency Room. Intensive Care Unit. Mobile Lab. Research Center

Marketing Department August 2022



Our History

OUR HISTORY

Murex started in 1998 by supplying top quality products to hospitals and clinics in Northern Luzon. At present, Murex is operating nationwide with dealers in Visayas and Mindanao and direct accounts in the National Capital Region. For the past two decades, Murex has become of the most trusted in vitro diagnostic products distributor in the Philippines by international and local partners.



2018

Murex celebrated its 20 years in service.

2022

Murex is currently in the ground works toward being ISO 9001:2015 certified





ON BUSINESS | To provide best in class collective products and medical devices for in vitro diagnostics



ON CUSTOMERS | To be present and accessible to our customers by establishing dealers in each region and ensuring customer delight in every interaction



ON BUSINESS PARTNERS | To sustain and achieve exclusive distributorship of our products and maintain firm, robust, and long-standing relationships with our business partners



ON EMPLOYEES | To create a training academy for continuous development, learning, and growth of the most valuable asset of our organisation – our people



ON SOCIETY | To promote literacy through continuous and sustainable education in the communities we are actively part of





OUR VISION

To be one of the **leading in vitro diagnostics medical devices distributors** nationwide by 2025

OUR CORPORATE VALUES

Value on Loyalty | We adhere to and abide by the cause and culture established by our MDPS founders

Value on Dedication & Hardwork | Our sincere and steadfast commitment to and passion for MDPS

Value on Teamwork | We transcend and translate our individual goals to the goals of MDPS through collaboration with others



Our Vision & Corporate Values



OUR QUALITY POLICY


Our commitment is to provide **quality products** with **excellent customer service** by our **loyal**, **dedicated**, **and hardworking team of people**, through our **efficient operational systems** and **effective communication lines**, to achieve **customer satisfaction** and **optimum organizational performance**





ax Clearar	ice ce	rtilicate
BUREAU OF	C OF THE PHILIPPINES RIMENT OF FIRMACE INTERNAL REVENT AL NO.1, CALASIAO, PANGASIN	
	TCBP NO. 01-006	05-04-R0487-2022-M
TAX CLEARA	NCE CERT nt to Executive Order No. 398)	
** JIMENEZ, M	ELISSA CE	NDANA **
* (MUREX DIAGNOSTI		
8 BONIFACIO ST., POBLAC	ION EAST, SAN NICOLAS	5, PANGASINAN
	628-398-00000 Identification Number	
This is to certify that the above mentione Certificate having satisfied all the otheria se to Revenue Regulations No. 8-2016, as am	t forth by the BIR as of the date	ence of this Tax Clearance e of this certification pursuant
Tax liabilities recorded after the aforesaid d by this tax dearance.	ates or outside the jurisdiction	of this Office are not covered
Issued this 4° day of May, 2022.		
NOTE THIS CERTIFICATE SHALL BE VALID AND IT INVOKED FOR VICKATION OF THE CITEDER SPEC INVOKED FOR VICKATION OF THE CITEDER SPEC INVOKET TRANSPORT OF THE VICKATION OF THE TRANSPORTION NO. BIE-COS22-DECOMPOSITION	I ED UNDER REVENUE REGULATION , WHICHEVEN COMES ZABUER. T COME FEE OF FLOR WAS PAID ON MA	S NO, 8 2006, AS AMEMORD AND INS NHALL NOT BE UNED ON IV 4, 2822 UNDER EFTS PAYMENT
	RIVIC G. TULIO	10580790
WARNING: Counterfeiting is punishat	Ne by Jaw. For outhestic	try, please with 818 website

DTI Registration



Business Permit



Philgeps Registration

Republic of the Philippine	3
Department of Budget and Man	agement
PROCUREMENT SERV	ICE
CERTIFICATE OF PHILGEPS	REGISTRATION
(Platinum Membershi	p)
THIS IS TO CERTIFY THAT	
MUREX DIAGNOSTIC PRODUCT	SSPECIALISTS
8 BONIFACIO STREET POBLACIO	
Son Nicolas , Pangasinan , Region I , I	Philippines
is registered in the Platlippine Government Electronic Procurem parsunt to Section 8.5.2 of the Revised Implementing Rules and observise known as the Government Procurement Reform Act.	
This further certifies that MUREX DIAGNOSTIC PRODUCTS SI eligibility documents in the PhilGEPS Supplier Registry as listed in Au and made an integral part heread.	ECIALISTS has submitted the require nex A, which document is associed here
For the purpose of apdating this Certificate, all Class 'A" eligibility Revised hyphemetring Rates and Regulations of Regulate, Act No. 980 validity of this Certificate which remain surverse and apdated? De fuilt Certificate with the carrier and updated Class 'A" eligibility docume of its validity and tasks that all of the experied Class 'A" eligibility	4 supporting the verticity, authenticity at re-by the prospective Bidder to update the stashall result in the automatic passessis
By submitting this Certificate, the Bidder certifies:	S
 the asthenticity, genaineness, validity, and completeness of th submitted; 	e copy of the original eligibility documen
2 the vertecity of the statements and information contained ther	ein:
 Heat the Certificate is not a guaranty that the named regists being determined to be such for that particular bidding, nor i the post-qualification stage; and 	
4. that any finding of concealment, falsification, or misreprese submitted, or the contents thereof shall be a ground for diraps bidding process, wildow prejudice to the imposition of app penalty in accordance with the lans.	slification from further participation in th
This Certificate is valid until 28-Apr-2028	
This is a system gen	Issued this 25th day of April 202 mated certificate. No signature is require
	(m) + 22 # 700 (R) 2 (R)
Documentary Stamp Tax Paid Php 30:00 Centificate Reference No: 200608107251883994731	Page 1 of 3



BIR 2303

License to Operate

R FORM 2303 EVISED APRIL 2010 R	EVENUE DIST	KAGAWARAN KAWANIHAN NG REVENUE REGION	NG PILIPINAS NG PANANALAPI RENYAS BITERNAS ING 001 - CALASIA 01- URDANETA, EA	0
				OCN: 006RC2022900060 Date OCN Generated: April 21, 3
	c	ERTIFICATE	OF REGISTI	RATION
TIN & BRANCH COD 908-428-390-0000 REGISTERING OFFIC REGISTERED ADDR	æ	NAME OF TAXPAY	A CENDANA	TIN ISSUANCE DATE August 8, 2007
		LAST 2447 BAN NIC	OLAS PANGASINAN	PHLIPPINES
TAX TYPES	FORM	FILING START DATE	FILING	FILING DUE DATE
	E 17010	February 21, 2013	QUARTERLY	1st Quarter-on or before MAY 2nd Quarter-on or before AUQUST 15 3rd Quarter-on or before November 15
INDIVIDUAL INCOM TAX	e 1701	Fabruary 21, 2013	ANNUALLY	On or before April 15 of each year covering income for the proceeding totalitie year.
REDISTRATION FE	E 0005	February 21, 2013	ANNUALLY	On or before the last day of January.
VALUE ADDED TA	x 2550M	August 22, 2014	MONTHLY	Not later than the 20th day following the close of the mon
WITHHOLDING TAX COMPENSATION		Fabruary 21, 2013	ANNUALLY	Dri schelsen January 11 of ih year following the colondar ye in which companisation payments and other income payments subject to grad withholding bas were paid or accurate.
WITHHOLDING TAX EXPANDED/OTHER		December 11. 2018	ANNUALLY	
WITHHOLDING TAX EXPANDED/OTHER	5 0010E	December 11, 2018	MONTHLY	On or before the 10th day of th month following the month in which with holding was made
WITHHOLDING TAX EXPANDED/OTHER		December 11. 2018	QUARTERLY	Not later than the later day of the month following the classe of the quarter during which withholds was made.
TAXPAYER TYPE/S	SNGLE	PHOPHETORISHIP	ONLY PRESIDENT O	
BUSINESS INFORM	TION DETAIL	s		
TRADE NAME 1	MUREX DIA	INDETIC PRODUCT	CATEO	IORY REGISTRATION DATE May 9, 2008
(P8IC)	47720-RETA PHARMACEI GOODS, CO		*	
Line of Business		OF OTHER GOOD		
(P\$#C)		R PERSONAL SERV	ICES Secon	dary

Republic of the Philippines Department of Health TODO AND DRUG ADMINISTRATION Concerning Concerning Concerning Administration Concerning
LICENSE TO OPERATE
Medical Device Importer BA
MUREX DIAGNOSTIC PRODUCTS SPECIALISTS #8 Bonifacio Street, Poblacion East, San Nicolas, Pangasinan
Overe:: MELISSA CENDAÑA IIMENEZ Licease Number: Application Type Variadan Validity of Licease: 27 December 2022
This LTO shall be reasoned within <u>90 days</u> before its expiration, upon compliance with relevant have, wices and regulations and the payment of fees. This LTO may be suspended, cancelled an enviced by this Offices for enause if found withining RA CPU fland reliable summites.
Eurofermone, with this LTO, the FDA allows the establishment to apply for a marker achiever and the end of the end of the stable and the stable probability from to manufacture, improvements and a confer for stable, distributions transfer, advertisement indice promition are the case may be
BY AUTHORITY OF THE DIRECTOR GENERAL
MARIA CECILIA C, MATIENZO Directo IV Center for Device Regulation Radiution Health and Research
A FBA
DISPLAY IN PUBLIC VIEW
BAT FBAT B FBAT B FBAT B FBAT B FBAT
FBA This electronic LTO is computer generated and their not require signature.
Verify at http://www.ide.on.ph.industry.comervedin-//70-tonion/restve2 FDA-0477831

Permits & Licenses



Bringing the latest IVD technology to the Philippines

Investing on travel and relationships, Murex Diagnostics was able to establish partnerships with manufacturers and suppliers from all over the world such as Europe, North America, and Asia.







MEDICAL DIAGNOSTICS EQUIPMENT

- CLINICAL CHEMISTRY
- HEMATOLOGY & COAGULATION
- IMMUNOLOGY
- CLINICAL MICROSCOPY
- OTHERS

VETERINARY DIAGNOSTICS EQUIPMEMT

- VETERINARY HEMATOLOGY
- VETERINARY BIOCHEMISTRY

OUR PRODUCTS







EMP-168

Semi-Automated Analyzer -----

- Windows based user friendly interface
- Support RS-232, USB, Ethernet, and SD card
- Light source: Halogen 6V/10W
- Repeatability: CV≤ 0.5%
- Carry Over $\leq 1\%$
- Linearity: r ≥ 0.999
- Flow cell: 32 µl with peltier temperature control (25, 30, 37°C)



STAT FAX 4500

Semi-Automated Analyzer

-
- Read well accepts 12 mm test tubes and 1 cm cuvettes
- Flow cell for rapid fluid sampling with low carryover
- 8 position filter wheel
- Standard 6 filter configuration (340 nm 700 nm)
- Optional read cell cooling at 25 ℃, 30 ℃, or 37 ℃
- Incubation block at 37℃
- Stores up to 120 tests in memory



Clinical Chemistry Analyzer (Wet Method)







ChemWell-TICS

Fully Automated Analyzer

Ease of Use

• ChemWell-T iCS software is intuitive. Visual prompts & graphics indicate placement of reagents/samples and assist the user with rapid setup. System includes a STAT facility, uses a multi constituent calibrator.

Liquid Stable Reagents

- No reagent prep work required. Load kit bottles directly onto
- yield long, stable calibrations.

Throughput

• 100 endpoint test per hour

EXC 200

Fully Automated Analyzer

User Friendly

- Colorful touch screen and intuitive user-friendly navigation menu
- Waste with high concentration and low concentration are discharged separately, more environmental friendly

Economic Usage

- Lower reaction volume: 90 µL
- Less water consumption: 5L/H
- Precise reagent absorption with step by 0.5 μL

Throughput

• 200T/hr



• the instrument. The reagent rack is removable for refrigeration; covered, cooled, and insulated while on board to

Clinical Chemistry Analyzer (Wet Method)





SAT 450

Fully Automated Analyzer

On board reagents New I.S.E. Module (OPTIONAL) **Cost Efficiency**

- cuvettes each) selective cuvettes wash with 3 different liquids:
- acid, alkaline and neutral;
- cuvettes quality check.

Throughput

- 280 photometric tests/hour + 160 ISE/hour
- up to 440 tests/hour (280 photometric tests/hour + 160 ISE/hour)



• 80 long life, 6 mm optical path, semi-disposable cuvettes (4 separate racks with 20

• low water consumption: 1.8 l/h; 200 µl minimum reaction volume; automatic

Clinical Chemistry Analyzer (Wet Method)



Reagent Name	Method	Packaging
Chema Albumin	Colorimetric EP	2 x 50 ml, 4x100 ml
Chema Alkaline Phos FL (DGKC)	DGKC, kinetic colorimetric	4x20 ml, 12x20 ml
Chema Alkaline Phos FL (DGKC)	IFCC, kinetic colorimetric	8x50 ml, 4 x 20 ml, 12x20 ml, 8x50 ml
Chema Amylase FL	CNP-G3, kinetic colorimetric	6x10 ml
Chema Bilirubin Total FL	DCA method	5x25 ml
Chema Bilirubin Direct FL	DCA method	5x25 ml
Chema Calcium ASX	Arsenazo (III) method	2x50 ml
Chema Chloride	Colorimetric EP	2x50 ml
Chema Cholesterol FL	Enzymatic colorimetric EP	2x50 ml, 4x100 ml, 4x250 ml



Clinical Chemistry Reagents (Wet Method)



Reagent Name	Method	Ра
Chema HDL - Direct FL	Immuno-inhibition	4x2
Chema HDL Prec Reagent	PEG precipitating reagent	4x´
Chema CK-MB FL DGKC/IFCC	IFCC/DGKC, kinetic UV	6x´
Chema Creatinine	Colorimetric fixed-time	4x´
Chema Creatinine-E FL	Enzymatic colorimetric EP	5x2
Chema Ethanol FL	Enzymatic UV	6x´
Chema Gamma GT FL	Szasz, kinetic colorimetric	4x2
Chema Glucose FL	Enzymatic colorimetric EP	4x´
Chema GOT-AST FL (IFCC)	IFCC, kinetic UV	4x2



ackaging

20 ml, 12x20 ml, 4x100 ml
(100 ml
:10ml
100ml, 4x125 ml, 4x250 ml
25 ml, 15x25 ml, 10x60 ml
(10 ml
20 ml
(100 ml, 4x250 ml
20 ml, 12x20 ml, 8x50 ml, 5x100 ml

Clinical Chemistry Reagents (Wet Method)



Reagent Name	Method	Ра
Chema GPT-ALT FL (IFCC)	IFCC, kinetic UV	4x2
Chema Iron CRX	Chromazurol B, colorimetric EP	2x5
Chema LDH FL (DGKC)	DGKC, kinetic UV	бх́
Chema LDL Direct FL	Immuno-inhibition	4 x
Chema Lipase FL	Kinetic colorimetric	бх́
Chema Magnesium XL blue	Xylidyl blue, colorimetric EP	6x!
Chema Phosphorus UV	UV endpoint	2x!
Chema Total Proteins	Biuret, colorimetric EP	2x!
Chema Triglycerides	Enzymatic colorimetric EP	2x!

Clinical Chemistry Reagents (Wet Method)

ackaging

20 ml, 12x20 ml, 8x50 ml, 5x100ml
:50 ml
(10 ml
x 20 ml
(10 ml
:50 ml
50 ml, 4x100ml
50 ml, 4x100ml
50 ml, 4x100 ml, 4x250 ml



Reagent Name	Method	Pa
Chema Urea UV FL	Fixed time UV	4x2
Chema Uric Acid T FL	Enzymatic colorimetric EP	4x
Chema Uric Acid AOX FL	Enzymatic colorimetric EP, AOX	5x5



ackaging

(20 ml, 12x20 ml, 8x50 ml, 5x100ml

x100 ml

<50 ml, 5x100ml

Clinical Chemistry Reagents (Wet Method)



Reagent Name	Method	Packaging
Chema Ethanol Control		2x5 ml
Chema Quantinorm normal Control	Human based control serum N	5x5 ml, 10x5 ml
Chema Quantipath normal Control	Human based control serum P	5x5 ml, 10x5 ml
Chema Multinorm	Human based control serum N	5x5 ml
Chema Multipath	Human based control serum P	5x5 ml
Chema Autocal H	Human based multicalibrator	5x3 ml, 10x3 ml

Clinical Chemistry Control & Calibrator (Wet Method)





Celercare CM5

Gen. Chem & Electroloyte Analyzer

- Xenon light source
- Microfluidic Reagent Rotor
- Sample-to-answer in 3 simple steps and approximately 7 minutes

Reagent Discs

	Reagent Name	Parameters
	Electrolyte Panel Lyophilized Kit	K+, Na+, Cl-, Ca2+, P, Mg2+, CO2
	General Chemistry I Lyophilized Kit	TP, ALB, GLO*, A/G*, ALT, AST, TBIL HDL-C, LDL-C*
	Clinical Emergency Lyophilized Kit	AST, UA, CRE, CK, CK-MB, LDH, a-HBD
	Liver and Renal Function Lyophilized Kit	TP, ALB, TBIL, ALT, AST, GGT, CRE, UR

.....

IL, DBIL, IBIL*, UA, CRE, UREA, GLU, TG, CHOL,

DH, AMY, GLU, K+, Na+, Cl-, CO2

JREA, GLU

Clinical Chemistry & Electrolyte Analyzer (Dry Method)





XI-931T

Electroloyte Analyzer

- High-performance electrodes: high precision and wide linearity range, with long shelf life (up to 12 months for Na, Cl electrodes).
- Support high/low-value calibrators and 3 levels of QC, automatic calibration of slope and intercept values, and display of whole QC graphs.
- Two sampling models: automatic sampling with an autosampler or manual sampling.
- Operational software can be upgraded through a USB port.
- Automatic liquid level detection and alarms. Real-time diagnostic of system status and reagent pack status information.

-	•	_	_		
1	-				C. Star
1	anne	101120		-	1.1.4.4
	1.11		24	1	

New generation of electrodes with 100% maintenance-free.

Auto sampler (option)

Na+	<
CI-	\leq
Ca2+	\leq
рН	\leq
CO2	<

Parameter

K+



CV%	Resolution	Range (mmol/L)
≤1.0%	0.01 mmol/L	0.50 - 15.00 mmol/L
≤ 1.0%	0.1 mmol/L	30.0 - 200.0 mmol/L
≤ 1.0%	0.1 mmol/L	30.0 - 200.0 mmol/L
≤ 3.0%	0.01 mmol/L	0.10 - 5.00 mmol/L
≤ 2.0%	0.01 mmol/L	4.0 - 9.50 (Unit)
≤ 3.0%	0.1 mmol/L	6.0 - 50.0 mmol/L

Electrolyte Analyzer (Wet Method)



Caretium Consumables & Control

Reagent Name	Packaging
Caretium Solution Pack	pack
Caretium K electrode	рс
Caretium Na electrode	рс
Caretium Cl electrode	рс
Caretium Ca electrode	рс
Caretium pH electrode	рс
Caretium Reference electrode	рс
Caretium internal electrode	рс
Caretium Ref internal electrode	рс
Caretium Reference column	рс



agı	ng

Electrolyte Analyzer (Wet Method)



Caretium Consumables & Control

Reagent Name	Packaging
Caretium QC solution Lv1	vial
Caretium QC solution Lv2	vial
Caretium QC solution Lv3	vial
Caretium Daily Cleaning Solution (100 ml)	bottle
Caretium Weekly Cleaning Solution (100 ml)	bottle
Caretium Na conditioner (100 ml)	bottle
Caretium Urine Solution	bottle



Electrolyte Analyzer (Wet Method)





GREENCARE A1c

One Step HbA1c measuring system

- Provide accurate quantitative test results for HbA1c
- NGSP and IFCC certified and CV is less than 3%
- Provide the most convenient solution
- Help clinician decision during the patient's visit
- Hemoglobin A1c Test Kit
 - Automatic recognition of code and expiration date
 - Store and operate at Room temperature (2-30 °C)

Reagent Name	
Greencare A1c Test Kit	107
Greencare A1c check kit	kit



ckaging

T/Box






DH26

3-Part Hematology Analyzer 21 Parameters Sample Volume: 9µL

- Compact design
- Reliable performance
- Humanized design
- Cost-effective innovation





DF51

5-Part Hematology Analyzer 23 Parameters Sample Volume: 20 ul

- Reliable component
- Valve from Japan
- Needle from Switzerland
- Liquid pump from Germany
- Long life syringe



DH73 3-Part & 5-Part Combined Auto Hematology Analyzer CBC+DIFF: 25 Parameters CBC+3DIFF: 21 Parameters Sample Volume: 20 ul

- 5-Part testing mode
- 3-Part testing mode
- Fully Automatic Intelligent Technology
- Excellent Data management system

Hematology Analyzer (Wet Method)



DH26 Reagent

Reagent Name	Packaging
DIL-K	20L
LYK-1	500ml
CLE-P	50ml

DH73 Reagent

Reagent Name	
DIL-A	
LYG-1	
LYA-2	
LYA-3	
CLE-P	

DF51 Reagent

Reagent Name	Packaging
DIL-C	20L
LYC-1	200ml
LYC-2	500ml
CLE-P	50ml

Control Blood

Reagent Name

DM-3D (H,M,L)

DM-5D (H,M,L)



Packaging	
20L	
500ml	
500ml	
1L	
50ml	

DH26, DH73

DF51, DH73

Hematology Analyzer (Wet Method)





Micropoint Q3

Coagulation

- Reliable laboratory results in 30-100 seconds for PT and 3-7 minutes for aPTT
- 10 ul of fingertip whole blood, fresh venous whole blood or citrate whole blood
- On board quality control
- No calibration or instrument manipulation required

Q3 Reagents & Control

Reagent Name
Qlabs PT-INR Test Strips (Pro)
Qlabs aPTT Test Strips
Qlabs Coag Panel 2 Test Strips
Qlabs PT-INR Liquid Control
Qlabs aPTT Liquid Control
Qlabs Coag Panel 2 Liquid Control
Qlabs Re-Calcification Reagent



.....

Packaging
24 strips
12 strips
12 strips
set
set
set
set





ichroma II

2 level calibration Up to 30 tests/ hour On-board QC & System check Battery and AC operated

- Built-in printer
- Built-in multi-timers





ichroma III

High throughput results Random access Multi cartridges testing Temperature control (25 – 35 °C) LIS/HIS connectivity

- Simple, Fast, Accurate
- Suitable for the clinical laboratories
- No need thermal chamber
- 12 tests available at once
- Up to 57 Test Parameters

Immunology Analyzer (FIA)



ichroma Test Parameters

Cardiac	Cancer	Hormone	Infection
Tn-I Tn-I Plus	PSA AFP	TSH TSH Plus	CRP PCT
CK-MB D-Dimer NT-proBNP	CEA iFOB Neo	T3 T4 FSH	ASO Dengue IgG Dengue Ns
Myoglobin hsCRP Cardiac Triple ST2	Diabetes HbA1c Microalbumin Cystatin C	Progesterone LH PRL Testosterone Cortisol AMH	COVID19 A COVID19 A IL6

n

Autoimmune

RF IgM Anti-CCP Plus

gG/lgM ls1 Ag Ab Ag

Others

Ferritin Vitamin D Calprotectin H. pylori SA


Immunology Analyzer (FIA)





iFlash 1200

..... Up to 500 Reaction vessels 16 sample positions 120 T/H 115 reagents position 120+ Broader assay menu **Effective differential diagnosis**





Immunology Analyzer (CLIA)



Reagent Name

QuickProfile Dengue IgG/IgM Combo Test Card

QuickProfile Dengue NS1 Antigen Test Card

QuickProfile Dengue IgG/IgM & NS1 Duo Test Card

QuickProfile S. typhi IgG/IgM Combo Test Card

QuickProfile Malaria pf/pan Antigen Test Card

QuickProfile H.pylori Stool Antigen Test Card

QuickProfile H.pylori Antibody Test Card

QuickProfile COVID-19 Antigen Test Strip

QuickProfile 2019 nCoV lgG/lgM Test Card



	•	
Pac	kaσı	nσ
Pac	N GSI	'' 8

25T/Box
25T/Box
10T/Box
25T/Box
25T/Box
25T/Box
40T/Box
20T/Box
25T/Box

Immunology (Rapid Test Division)



Reagent Name

QuickProfile HIV I&II Test Card*

QuickProfile HBsAg Test Card*

QuickProfile HAV IgM Test Card

QuickProfile Syphilis Antibody Test Card*

QuickProfile MET/ THC Test Card (MET1000/THC500)*

QuickProfile Cardiac Panel Test Card (Tnl/CK-MB/Myoglobin)

QuickProfile Troponin I Test Card

QuickProfile Fecal Occult Blood Test Card



Packaging

40T/Box
40T/Box
25T/Box
40T/Box
40T/Box
25T/Box
40T/Box
25T/Box

Immunology (Rapid Test Division)



Reagent Name

Basechek Multi-Drug 5 Drugs Rapid Test Cassette (Urine)

Basechek Leptospira IgG/IgM Rapid Test Cassette

Basechek Multi-Drug 2 Drugs Rapid Test Cassette (Urine) MET 500/ THC50

Basechek hCG Pregnancy Urine Test Cassette



—	•
Packa	aging
	69

25T/Box

25T/Box

25T/Box

25T/Box

Immunology (Rapid Test Division)





Uryxxon Relax

Up to 500 Reaction vessels 16 sample positions 120 T/H 115 reagents position 120+ Broader assay menu Effective differential diagnosis



Reagent Name

Uryxxon Stick 10

Medi-Test Contr



Uryxxon 500

Quick - Reads 400 strips per hour Easy - Touch screen operation Hygienic - Automatic strip detection **Sediment**

- 30 user defined sediment parameters with 5 levels each
- Optional external sediment keyboard



Uryxxon Stick 10

(Glucose, Ketones, Protein/ Albumin/ Blood/ Leukozytes, Nitrile, Density/ pH value/ Bilirubin/ Urobiliongen

e	Packaging
0	100
rol	set



Metro Manila Office:

119 14th Ave., Cubao 1109, Quezon City, Metro Manila T/F +02-8-962-8743

26 East Road, San Roque, Quezon City, Metro Manila

Pangasinan Office:

8 Bonifacio St., San Nicolas 2447, Pangasinan T/F +075-600-8159 +075-632-4435

THANK YOU

info@murexdiagnostics.com **f** @murexdiagnostics.ph